



Testing Laboratory Basics

Below are some answers to common questions we receive for testing at the BioTracking testing lab. If the information listed below does not answer your questions, additional information can be found by visiting www.biotracking.com and clicking on BioTracking Testing Laboratory Information. If you have any additional questions, please contact us by email at testinglab@biotracking.com or phone at (208) 882-9736. Thank you!

Test Name	Sample Volume & Collection Vial	Coolant During Shipping	Recommended Shipping Method	Testing Requirements	Price per Sample
BioPRYN[®] Cattle ELISA	2 ccs of blood in a red top tube or 1 cc of serum	No coolant is required	USPS, FedEx or UPS	28 (cows) or 25 (heifers) days post breeding (by AI or bull removal) or 25 days post implant (32 days age of embryo) & 73 days post calving	\$2.50
BioPRYN[®] QK Cattle ELISA	2 ccs of blood in a red top tube or 1 cc of serum	No coolant is required	USPS, FedEx or UPS	28 (cows) or 25 (heifers) days post breeding (by AI or bull removal) or 25 days post implant (32 days age of embryo) & 90 days post calving	\$3.50
BioPRYN[®] Sheep & Goat ELISA	2 ccs of blood in a red top tube or 1 cc of serum	No coolant is required	USPS, FedEx or UPS	30 days post breeding (by AI or buck/ram removal)	\$6.50
BioPRYN[®] Bison ELISA	2 ccs of blood in a red top tube or 1 cc of serum	No coolant is required	USPS, FedEx or UPS	40 days post breeding (by AI or bull removal)	\$5.00
BioPRYN[®] wild Wildlife ELISA	2 ccs of blood in a red top tube or 1 cc of serum	No coolant is required	USPS, FedEx or UPS	40 days post breeding (by AI or buck/bull removal)	\$20.00
BioPRYN[®] es Equine ELISA	2 ccs of blood in a red top tube or 1 cc of serum	Coolant is required	FedEx or UPS Overnight delivery (blood) 2 nd day delivery (serum)	70 days post breeding (by AI or stud removal) & prior to 2 weeks of foaling	\$25.00
BVD ELISA	2 ccs of blood in a red top tube or 1 cc of serum or 1 cm by 1cm ear notch tissue	Coolant is required during summer months for both ear notches and blood	USPS, FedEx or UPS (Samples should not be in transit for more than a couple days)	Blood or serum sample: prior to colostrum or 3 months of age or older Ear notch sample: any age	\$3.65
CAE/OPP ELISA	2 ccs of blood in a red top tube or 1 cc of serum	Coolant is required during summer months	USPS, FedEx or UPS (Samples should be tested within 2 weeks of blood draw)	9 months of age or older & 6-8 weeks after vaccinations/illness or injury	\$4.00